

Sodium Test Kit

REF	Pack Size	Reagent 1	Reagent 2 (Std)
SODLMS01	25x1 ml	25x1 ml	1x2 ml

INTENDED USE

Sodium reagent is used for quantitative estimation of Sodium concentration in human serum or plasma

CLINICAL SIGNIFICANCE

Sodium is the major positive ion in fluids outside of cells. Most of the sodium in the body (about 85%) is found in blood and lymph fluid. It ensures a proper fluid and electrolyte or pH balance in our body, together with chloride and potassium. It enables our cell wall to draw in nutrients. It plays a role in nerve function and muscle contraction, in controlling the heartbeat by helping in its origin and maintenance. Sodium levels in the body are the adrenal glands. Aldosterone levels tell the kidney when to hold sodium in the body instead of passing it in the urine. Too much or too little sodium therefore can cause cells to malfunction, and extremes in the blood sodium levels (too much or too little) can be fatal. Too much sodium in the diet may raise blood pressure in some people. For those who have high blood pressure, eating foods with a lot of sodium makes their chance of heart disease, stroke, and kidney damage higher, heart failure gets worse when too much sodium is eaten. It increases the amount of water the body holds in and this causes swelling of the legs and hands. Increased sodium (hypernatremia) in the blood occurs wherever there is excess sodium in relation to water. There are numerous causes of hypernatremia: these may include kidney disease, too little water intake, and loss of water due to diarrhea and/or vomiting. High levels of sodium in the body are associated with high blood pressure and hypertension. Low sodium levels are uncommon and most often occur as a side effect of taking medicines that make you urinate more, such as diuretics. Severe diarrhea or vomiting or heavy sweating may also cause low sodium levels.

METHOD: Colorimetric method

TEST PRINCIPLE

The method is based on reaction of sodium with a selective chromogen producing a chromophore whose absorbance varies directly as the concentration of sodium in the test specimen.

KIT CONTENTS AND COMPONENTS

Reagent 1 : Sodium Reagent
Reagent 2 : Standard 150mmol/L

MATERIAL REQUIRED BUT NOT PROVIDED

Laboratory instrumentation, spectrophotometer UV/VIS with thermostatic cuvette holder or clinical chemistry analyzer: semi-automated, calibrated micropipettes, glass or high-quality polystyrene cuvettes, test tubes/rack, heating bath, controls, saline.

SAFETY PRECAUTIONS AND WARNING

- For in-vitro diagnostics use only.
- Do not pipette by mouth. Avoid contact with skin and eyes. If spilt thoroughly wash affected area with water.
- Do not use the reagent after the expiration date printed on the kit.

REAGENT PREPARATION STORAGE AND STABILITY

The reagent is ready to use.
The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C.. Protect from direct sunlight.

REAGENT DETERIORATION

Discard the reagent if absorbance exceeds 1.0 against distilled water.

PROGRAM

Reaction Mode	End point
Wavelength	630 nm
Light Path	10 mm
Blanking	Reagent blank
Reagent Volume	1000 µl
Standard Volume	10 µl

Sample Volume	10 µl
Incubation	10 mins at 37°C
Standard Conc.	150 mmol/L
Linearity	180 mmol/L

PROCEDURE

Reagent	Blank	Standard	Test
R1 Sodium Reagent	1000 µl	1000 µl	1000 µl
Distilled water	10 µl		
R2- Standard		10 µl	
Sample			10 µl

Mix well and incubate for 10 minutes at 37°C. Read the absorbance of test sample A(T) and standard sample A(S) against reagent blank (B)

CALCULATION

Sodium concentration = $A(T) / A(S) \times 150 \text{ mmol/L}$ (standard concentration)

NORMAL VALUES

135-155 mmol/L

It is recommended for each laboratory to establish its own reference ranges for local population.

LIMITATIONS

If the value exceeds 180mmol/L dilute the sample with 0.9% saline solution and result multiplied by dilution factor.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform internal quality control with assayed normal and assayed abnormal, to confirm the validity of the test and assure the accuracy of patient result. Using the recommended calibrator or the Standard included, calibrate the assay.

WASTE DISPOSAL










This Product is made to be used in professional laboratories.

HIGHLIGHTS

- Storage condition mentioned on the kit is to be used.
- Before testing bring the reagents to the R.T.
- Avoid reagents contamination.
- Every time use new pipette-tips for pipetting out the reagents.
- These Reagent kits meant for laboratory/professional use only, not for Drug use.

REFERENCE

Tietz N.W., ed. Clinical Guide to laboratory Tests, 3rd ed. Philadelphia, Pa: W.B. Saunders, 1995:466;487
Lords Data File.

	Catalog No.		Contain Sufficient for test
	Batch No.		Instruction for use
	Manufacturing Date		In-vitro Diagnostics
	Expiry Date		Storage temperature
	Manufacturer		

IFU/SOD/01 Rev.: 02; Rev Dated.: 22/07/2024

